


BTPE-10

PE-10 tubing, .011x.024in, 30m spool, non-sterile



SPECIFICATIONS

Tube material	Low density polyethylene
Inner diameter	.011±.002in (.23-.33mm)
Outer diameter	.024±.002in (.56-.66mm)
Nominal volume	0.6µL/cm (can range from 0.4-0.8 with ID variability)
Length on spool	30m (98 +5/-0ft)
Durometer	45D
Spool dimensions	.9x4.8in (2.2x12.2cm)
Sterility	Non-sterile (sterilize by ETO; do not autoclave)

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/tubing-connectors-pinports/tubing/polyethylene>

Instructions For Use

Product Video

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com