


BTPE-60

PE-60 tubing, fits 20ga, .030x.048in, 30m spool, non-sterile



SPECIFICATIONS

Tube material	Low density polyethylene
Inner diameter	.030±.002in (.71-.81mm)
Outer diameter	.048±.002in (1.2-1.3mm)
Nominal volume	4.5µL/cm (can vary from 3.9-5.1 with ID variability)
Length on spool	30m (98 +5/-0ft)
Durometer	45D
Spool dimensions	.9x4.8in (2.2x12.2cm)
Sterility	Non-sterile (sterilize by ETO; do not autoclave)

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/tubing-connectors-pinports/tubing/polyethylene>

Instructions For Use

Product Video

NOTE: All products are designed and sold for research only; not for human use.
