


BTPE-90

PE-90 tubing, .034x.050in, 30m spool, non-sterile



SPECIFICATIONS

Tube material	Low density polyethylene
Inner diameter	.034±.002in (.81-.91mm)
Outer diameter	.050±.002in (1.21-1.32mm)
Nominal volume	5.8μL/cm (can vary from 5.1-6.5 with ID variability)
Length on spool	30m (98 +5/-0ft)
Durometer	45D
Spool dimensions	.9x4.8in (2.2x12.2cm)
Sterility	Non-sterile (sterilize by ETO; do not autoclave)

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/tubing-connectors-pinports/tubing/polyethylene>

Instructions For Use

Product Video

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com