


Part Number: **C19PU-HFV1706**



not to scale

Full Description: Catheter for hamster femoral vein, PU 1.9Fr 4.5cm, 3Fr 25cm, collars @ 3.0, 4.5, 7.5cm. Fits 22ga.

Unit: each **Pouch Desc:** Hamster FVC PU 1.9-3Fr 30cm **Pouch:** 4x10

Approx Dead Vol: 80 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: hamster **Vessel:** femoral vein

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 4.5 cm **L2:** 25 cm
Material (IDxOD): PU 1.9Fr (.31x.64 \pm .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 \pm .03mm)

FEATURES

| Distal (tip) end | | Location | Fixed |
|------------------|--|-----------|-------------------------------------|
| Feature A: | Small silicone collar for 2Fr (0.9mm OD), 1mm wide | a: 3.0 cm | <input type="checkbox"/> |
| Feature B: | Large silicone collar for 3Fr (1.7mm OD), 3mm wide | b: 4.5 cm | <input checked="" type="checkbox"/> |
| Feature C: | Small silicone collar for 3Fr (1.3mm OD), 2mm wide | c: 7.5 cm | <input type="checkbox"/> |

| Proximal end | | Location | Fixed |
|--------------|--|----------|--------------------------|
| Feature Z: | | z: cm | <input type="checkbox"/> |
| Feature Y: | | y: cm | <input type="checkbox"/> |
| Feature X: | | x: cm | <input type="checkbox"/> |

PARTS ADDED TO POUCH

Part 1: **Part 2:**