


Part Number: C20PU-MJV2013



not to scale

Full Description: Mouse JVC, PU 2Fr 10cm, 13mm insertion length (a1.2cm), round tip. LS25 and SP25/12. Fits 25ga.

Unit: each

Pouch Desc: Mouse JVC 2Fr PU L10cm 13mm

Pouch: 3x7.5

Approx Dead Vol: 15 μ L

Connects to: 25ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 10 cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

	Location	Fixed
Distal (tip) end		
Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 1.2 cm	☒
Feature B:	b: cm	☐
Feature C:	c: cm	☐
Proximal end		
Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1: LS25

Part 2: SP25/12