

VABM1B/25

Vascular Access Button for mouse, magnetic, 1 channel 25ga, injector



SPECIFICATIONS

Number of ports	1
Port dead volume	3µL
Catheter connector	25ga (.020in, .51mm OD)
Compatible catheters	2Fr PU with ID of .41-.43mm
Syringe access	With PNP3M injector only (do not use sharp or blunt needles)
Materials in fluid path	Acetal, stainless steel, silicone
Compatible tethers	Instech VABM1T/25 or KVABM1T/25
Port pressure rating	150PSI (do not exceed 4lbs/2kg on plunger of 1ml syringe)
Felt material	Surgical grade polyester
Felt diameter	.56in (14mm)
Button weight	0.3g
Sterilization	ETO gas, one time use
Intended use	Exteriorization of mouse catheters

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/buttons-tethers/vascular-access/mouse-one-channel>

Instructions For Use <https://www.instechlabs.com/hubfs/pdfs/manuals/ifu-vabm.pdf>

Product Video <https://www.instechlabs.com/view-video-mouse-vab-use>

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com