

VABR1B/22

Vascular Access Button for rat, magnetic, 1 channel 22ga (white port), injector



SPECIFICATIONS

Number of ports	1
Port dead volume	4μL
Catheter connector	22ga (.028in, .71mm OD)
Compatible catheters	3Fr PU with ID of .60-.64mm
Syringe access	With PNP3M injector only (do not use sharp or blunt needles)
Materials in fluid path	Acetal, stainless steel, silicone
Compatible tethers	Instech VABR1T or VABR2T tethers
Port pressure rating	150PSI (do not exceed 4lbs/2kg force on plunger of 1ml syringe)
Felt material	Surgical grade polyester
Felt diameter	1in (25mm)
Button weight	1.9g
Sterilization	ETO gas, one time use
Intended use	Exteriorization of rat catheters

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/buttons-tethers/vascular-access/rat-one-channel>

Instructions For Use <https://www.instechlabs.com/hubfs/pdfs/manuals/VABR.pdf>

Product Video <https://www.instechlabs.com/view-video-rat-vab-use>

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com